

COLQ Antibody

Catalog No: #45327

Package Size: #45327-1 50ul #45327-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	COLQ Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	COLQ Antibody detects endogenous levels of total COLQ
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human COLQ
Target Name	COLQ
Other Names	COLQ, AChE Q subunit, EAD
Accession No.	Swiss-Prot#: Q9Y215NCBI Gene ID:8292
Uniprot	Q9Y215
GeneID	8292;
Calculated MW	47 kDa
Concentration	1.0mg/mL
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis of COLQ expression in SP20 (A), H9C2 (B), PC3 (C) whole cell lysates.

Note: This product is for in vitro research use only