

Defensin β 1 Antibody

Catalog No: #45331

Package Size: #45331-1 50ul #45331-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Defensin β 1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	Defensin β 1 Antibody detects endogenous levels of total Defensin β 1
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human Defensin β 1
Target Name	Defensin β 1
Other Names	DEFB1, BD1, Beta-defensin 1, BD-1, Beta Defensin 1, DEFB-1, DEFB101, Defensin, beta 1, Defensin beta, HBD-1, Beta-defensin-1, HBD1
Accession No.	Swiss-Prot#: P60022NCBI Gene ID:1672
Uniprot	P60022
GeneID	1672;
Calculated MW	7 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis of Defensin beta 1 expression in HEK293T (A), mouse liver (B), rat heart (C) whole cell lysates.

Note: This product is for in vitro research use only