

DIO3 Antibody

Catalog No: #45333

Package Size: #45333-1 50ul #45333-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	DIO3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	DIO3 Antibody detects endogenous levels of total DIO3
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human DIO3
Target Name	DIO3
Other Names	Type III iodothyronine deiodinase; 5DIII; DIOIII; Type 3 DI; Type-III 5'-deiodinase; DIO3; ITDI3; TXDI3
Accession No.	Swiss-Prot#: P55073NCBI Gene ID:1735
Uniprot	P55073
GeneID	1735;
Calculated MW	31 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details


Western blotting: 1:1000-3000

Images



Western blot analysis of lysate from HUVEC cells, using DIO3 antibody.

Western Blot analysis of various cells using DIO3 Polyclonal Antibody



Note: This product is for in vitro research use only