

## EAAT2 Antibody

Catalog No: #45336



Package Size: #45336-1 50ul #45336-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	EAAT2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	EAAT2 Antibody detects endogenous levels of total EAAT2
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human EAAT2
Target Name	EAAT2
Other Names	SLC1A2, EAAT-2, GLT-1, Glutamate transporter 1, EAAT2, GLT1
Accession No.	Swiss-Prot#: P43004NCBI Gene ID:6506
Uniprot	P43004
GeneID	6506;
Calculated MW	62 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000-3000

Note: This product is for in vitro research use only