

Inhibin β -C Antibody

Catalog No: #45359

Package Size: #45359-1 50ul #45359-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Inhibin β -C Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	Inhibin β -C Antibody detects endogenous levels of total Inhibin β -C
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human Inhibin β -C
Target Name	Inhibin β -C
Other Names	INHBC, Activin beta-C chain, IHBC, Inhibin, beta C, Inhibin beta c chain
Accession No.	Swiss-Prot#: P55103NCBI Gene ID:3626
Uniprot	P55103
GeneID	3626;
Calculated MW	38 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis of hela lysis using INHBC antibody.
Secondary antibody was diluted at 1:20000

Note: This product is for in vitro research use only