

Cleaved-Caspase-9 (D353) Polyclonal Conjugated Antibody

Catalog No: #C40503

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Package Size: #C40503-AF350 100ul #C40503-AF405 100ul #C40503-AF488 100ul

#C40503-AF555 100ul #C40503-AF594 100ul #C40503-AF647 100ul

#C40503-AF680 100ul #C40503-AF750 100ul #C40503-Biotin 100ul

Description

Product Name	Cleaved-Caspase-9 (D353) Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Specificity	Cleaved-Caspase-9 (D353) Polyclonal Antibody detects endogenous levels of fragment of activated Caspase-9 protein resulting from cleavage adjacent to D353.
Immunogen Description	Synthesized peptide derived from the C-terminal region of human Caspase-9.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CASP9; MCH6; Caspase-9; CASP-9; Apoptotic protease Mch-6; Apoptotic protease-activating factor 3; APAF-3; ICE-like apoptotic protease 6; ICE-LAP6
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	17
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only