

## CRF-RI Polyclonal Conjugated Antibody

Catalog No: #C40785



Package Size: #C40785-AF350 100ul #C40785-AF405 100ul #C40785-AF488 100ul  
 #C40785-AF555 100ul #C40785-AF594 100ul #C40785-AF647 100ul  
 #C40785-AF680 100ul #C40785-AF750 100ul #C40785-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
 Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CRF-RI Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Specificity	CRF-RI Polyclonal Antibody detects endogenous levels of CRF-RI protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human CRF-RI.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CRHR1; CRFR; CRFR1; CRHR; Corticotropin-releasing factor receptor 1; CRF-R-1; CRF-R1; CRFR-1; Corticotropin-releasing hormone receptor 1; CRH-R-1; CRH-R1
Accession No.	Swiss-Prot#:P34998NCBI Gene ID:1394
Uniprot	P34998
GeneID	104909134;1394;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	51
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

## Application Details

## Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

---

Note: This product is for in vitro research use only