

CRF-RII Polyclonal Conjugated Antibody

Catalog No: #C40786



Package Size: #C40786-AF350 100ul #C40786-AF405 100ul #C40786-AF488 100ul
 #C40786-AF555 100ul #C40786-AF594 100ul #C40786-AF647 100ul
 #C40786-AF680 100ul #C40786-AF750 100ul #C40786-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	CRF-RII Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Specificity	CRF-RII Polyclonal Antibody detects endogenous levels of CRF-RII protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human CRF-RII.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CRHR2; CRF2R; CRH2R; Corticotropin-releasing factor receptor 2; CRF-R-2; CRF-R2; CRFR-2; Corticotropin-releasing hormone receptor 2; CRH-R-2; CRH-R2
Accession No.	Swiss-Prot#:Q13324NCBI Gene ID:1395
Uniprot	Q13324
GeneID	1395;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	47
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only