

Dyrk1A Polyclonal Conjugated Antibody

Catalog No: #C40858



Package Size: #C40858-AF350 100ul #C40858-AF405 100ul #C40858-AF488 100ul
 #C40858-AF555 100ul #C40858-AF594 100ul #C40858-AF647 100ul
 #C40858-AF680 100ul #C40858-AF750 100ul #C40858-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------------|--|
| Product Name | Dyrk1A Polyclonal Conjugated Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Specificity | Dyrk1A Polyclonal Antibody detects endogenous levels of Dyrk1A protein. |
| Immunogen Description | Synthesized peptide derived from the N-terminal region of human Dyrk1A. |
| Conjugates | Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750 |
| Other Names | DYRK1A; DYRK; MNB; MNBH; Dual specificity tyrosine-phosphorylation-regulated kinase 1A; Dual specificity YAK1-related kinase; HP86; Protein kinase minibrain homolog; MNBH; hMNB |
| Accession No. | Swiss-Prot#:Q13627NCBI Gene ID:1859 |
| Uniprot | Q13627 |
| GeneID | 1859; |
| Excitation Emission | AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm |
| Calculated MW | 90 |
| Formulation | 0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide |
| Storage | Store at 4°C in dark for 6 months |

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only