

## Elf-1 Polyclonal Conjugated Antibody

Catalog No: #C40874



Package Size: #C40874-AF350 100ul #C40874-AF405 100ul #C40874-AF488 100ul  
 #C40874-AF555 100ul #C40874-AF594 100ul #C40874-AF647 100ul  
 #C40874-AF680 100ul #C40874-AF750 100ul #C40874-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
 Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Elf-1 Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Specificity	Elf-1 Polyclonal Antibody detects endogenous levels of Elf-1 protein.
Immunogen Description	Synthesized peptide derived from the C-terminal region of human Elf-1.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	ELF1; ETS-related transcription factor Elf-1; E74-like factor 1
Accession No.	Swiss-Prot#:P32519NCBI Gene ID:1997
Uniprot	P32519
GeneID	1997;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	67
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

## Application Details

## Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

---

Note: This product is for in vitro research use only