

## GPR56 Polyclonal Conjugated Antibody

Catalog No: #C40984



Package Size: #C40984-AF350 100ul #C40984-AF405 100ul #C40984-AF488 100ul  
 #C40984-AF555 100ul #C40984-AF594 100ul #C40984-AF647 100ul  
 #C40984-AF680 100ul #C40984-AF750 100ul #C40984-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
 Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	GPR56 Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Specificity	GPR56 Polyclonal Antibody detects endogenous levels of GPR56 protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human GPR56.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	GPR56; TM7LN4; TM7XN1; G-protein coupled receptor 56; Protein TM7XN1
Accession No.	Swiss-Prot#:Q9Y653NCBI Gene ID:9289
Uniprot	Q9Y653
GeneID	9289;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	78
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

## Application Details

## Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

---

Note: This product is for in vitro research use only