

Histone 1.0 Polyclonal Conjugated Antibody

Catalog No: #C41008



Package Size: #C41008-AF350 100ul #C41008-AF405 100ul #C41008-AF488 100ul
 #C41008-AF555 100ul #C41008-AF594 100ul #C41008-AF647 100ul
 #C41008-AF680 100ul #C41008-AF750 100ul #C41008-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Histone 1.0 Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Specificity	Histone 1.0 Polyclonal Antibody detects endogenous levels of Histone 1.0 protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human Histone 1.0.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	H1F0; H1FV; Histone H1.0; Histone H1'; Histone H1(0)
Accession No.	Swiss-Prot#:P07305NCBI Gene ID:3005
Uniprot	P07305
GeneID	3005;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	20
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only