

Histone deacetylase 10 Polyclonal Conjugated Antibody

Catalog No: #C41010

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #C41010-AF350 100ul #C41010-AF405 100ul #C41010-AF488 100ul

#C41010-AF555 100ul #C41010-AF594 100ul #C41010-AF647 100ul

#C41010-AF680 100ul #C41010-AF750 100ul #C41010-Biotin 100ul

Description

| | |
|-----------------------|--|
| Product Name | Histone deacetylase 10 Polyclonal Conjugated Antibody |
| Host Species | Rabbit |
| Clonality | Polyclonal |
| Specificity | Histone deacetylase 10 Polyclonal Antibody detects endogenous levels of Histone deacetylase 10 protein. |
| Immunogen Description | Synthesized peptide derived from the N-terminal region of human Histone deacetylase 10. |
| Conjugates | Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750 |
| Other Names | HDAC10; Histone deacetylase 10; HD10 |
| Accession No. | Swiss-Prot#:Q969S8NCBI Gene ID:83933 |
| Uniprot | Q969S8 |
| GeneID | 83933; |
| Excitation Emission | AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm |
| Calculated MW | 75 |
| Formulation | 0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide |
| Storage | Store at 4°C in dark for 6 months |

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only