

iPLA2 γ Polyclonal Conjugated Antibody

Catalog No: #C41072



Package Size: #C41072-AF350 100ul #C41072-AF405 100ul #C41072-AF488 100ul
 #C41072-AF555 100ul #C41072-AF594 100ul #C41072-AF647 100ul
 #C41072-AF680 100ul #C41072-AF750 100ul #C41072-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	iPLA2 γ Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Specificity	iPLA2 γ Polyclonal Antibody detects endogenous levels of iPLA2 γ protein.
Immunogen Description	Synthesized peptide derived from the C-terminal region of human iPLA2 γ .
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	PNPLA8; IPLA22; IPLA2G; BM-043; Calcium-independent phospholipase A2-gamma; Intracellular membrane-associated calcium-independent phospholipase A2 gamma; iPLA2-gamma; PNPLA-gamma; Patatin-like phospholipase domain-containing protein 8; iPLA
Accession No.	Swiss-Prot#:Q9NP80NCBI Gene ID:50640
Uniprot	Q9NP80
GeneID	50640;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	88
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250
 AF405 conjugated: most applications: 1: 50 - 1: 250
 AF488 conjugated: most applications: 1: 50 - 1: 250
 AF555 conjugated: most applications: 1: 50 - 1: 250
 AF594 conjugated: most applications: 1: 50 - 1: 250
 AF647 conjugated: most applications: 1: 50 - 1: 250
 AF680 conjugated: most applications: 1: 50 - 1: 250
 AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only