

Lfc Polyclonal Conjugated Antibody

Catalog No: #C41105



Package Size: #C41105-AF350 100ul #C41105-AF405 100ul #C41105-AF488 100ul
 #C41105-AF555 100ul #C41105-AF594 100ul #C41105-AF647 100ul
 #C41105-AF680 100ul #C41105-AF750 100ul #C41105-Biotin 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Lfc Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Specificity	Lfc Polyclonal Antibody detects endogenous levels of Lfc protein.
Immunogen Description	Synthesized peptide derived from human Lfc around the non-phosphorylation site of S885.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	ARHGEF2; KIAA0651; LFP40; Rho guanine nucleotide exchange factor 2; Guanine nucleotide exchange factor H1; GEF-H1; Microtubule-regulated Rho-GEF; Proliferating cell nucleolar antigen p40
Accession No.	Swiss-Prot#:Q92974NCBI Gene ID:9181
Uniprot	Q92974
GeneID	9181;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	111
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only