

L-type Ca⁺⁺ CP γ 7 Polyclonal Conjugated Antibody

Catalog No: #C41113



Package Size: #C41113-AF350 100ul #C41113-AF405 100ul #C41113-AF488 100ul
 #C41113-AF555 100ul #C41113-AF594 100ul #C41113-AF647 100ul
 #C41113-AF680 100ul #C41113-AF750 100ul #C41113-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	L-type Ca ⁺⁺ CP γ 7 Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Specificity	L-type Ca ⁺⁺ CP γ 7 Polyclonal Antibody detects endogenous levels of L-type Ca ⁺⁺ CP γ 7 protein.
Immunogen Description	Synthesized peptide derived from the C-terminal region of human L-type Ca ⁺⁺ CP γ 7.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CACNG7; Voltage-dependent calcium channel gamma-7 subunit; Neuronal voltage-gated calcium channel gamma-7 subunit; Transmembrane AMPAR regulatory protein gamma-7; TARP gamma-7
Accession No.	Swiss-Prot#:P62955NCBI Gene ID:59284
Uniprot	P62955
GeneID	59284;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	40
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only