

## MDMX Polyclonal Conjugated Antibody

Catalog No: #C41128



Package Size: #C41128-AF350 100ul #C41128-AF405 100ul #C41128-AF488 100ul  
 #C41128-AF555 100ul #C41128-AF594 100ul #C41128-AF647 100ul  
 #C41128-AF680 100ul #C41128-AF750 100ul #C41128-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
 Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	MDMX Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Specificity	MDMX Polyclonal Antibody detects endogenous levels of MDMX protein.
Immunogen Description	Synthesized peptide derived from human MDMX around the non-phosphorylation site of S367.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	MDM4; MDMX; Protein Mdm4; Double minute 4 protein; Mdm2-like p53-binding protein; Protein Mdmx; p53-binding protein Mdm4
Accession No.	Swiss-Prot#:O15151NCBI Gene ID:4194
Uniprot	O15151
GeneID	4194;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	45
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

## Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

---

Note: This product is for in vitro research use only