

Na<sup>+</sup> CP type IX $\alpha$  Polyclonal Conjugated Antibody

Catalog No: #C41200



Package Size: #C41200-AF350 100ul #C41200-AF405 100ul #C41200-AF488 100ul  
 #C41200-AF555 100ul #C41200-AF594 100ul #C41200-AF647 100ul  
 #C41200-AF680 100ul #C41200-AF750 100ul #C41200-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
 Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Na <sup>+</sup> CP type IX $\alpha$ Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Specificity	Na <sup>+</sup> CP type IX $\alpha$ Polyclonal Antibody detects endogenous levels of Na <sup>+</sup> CP type IX $\alpha$ protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human Na <sup>+</sup> CP type IX $\alpha$ .
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	SCN9A; NENA; Sodium channel protein type 9 subunit alpha; Neuroendocrine sodium channel; hNE-Na; Peripheral sodium channel 1; PN1; Sodium channel protein type IX subunit alpha; Voltage-gated sodium channel subunit alpha Nav1.7
Accession No.	Swiss-Prot#:Q15858NCBI Gene ID:6335
Uniprot	Q15858
GeneID	6335;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	220
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

## Application Details

## Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250  
 AF405 conjugated: most applications: 1: 50 - 1: 250  
 AF488 conjugated: most applications: 1: 50 - 1: 250  
 AF555 conjugated: most applications: 1: 50 - 1: 250  
 AF594 conjugated: most applications: 1: 50 - 1: 250  
 AF647 conjugated: most applications: 1: 50 - 1: 250  
 AF680 conjugated: most applications: 1: 50 - 1: 250  
 AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

---

---

Note: This product is for in vitro research use only