

Tensin-1 Antibody

Catalog No: #45400

Package Size: #45400-1 50ul #45400-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Tensin-1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	Tensin-1 Antibody detects endogenous levels of total Tensin-1
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human Tensin-1
Target Name	Tensin-1
Other Names	TNS1, MST091, MST127, MSTP127, MST122, Tensin 1, MSTP091, Tensin, Tensin-1, MSTP122, MXRA6, TNS
Accession No.	Swiss-Prot#: Q9HBL0NCBI Gene ID:7145
Uniprot	Q9HBL0
GeneID	7145;
Calculated MW	185 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C


Application Details

Western blotting: 1:1000-3000

Images



Western Blot analysis of K562 cells using Tensin-1 Polyclonal Antibody



Western blot analysis of lysate from K562 cells, using Tensin-1 antibody.

Note: This product is for in vitro research use only