

## NK-1R Polyclonal Conjugated Antibody

Catalog No: #C41243



Package Size: #C41243-AF350 100ul #C41243-AF405 100ul #C41243-AF488 100ul  
 #C41243-AF555 100ul #C41243-AF594 100ul #C41243-AF647 100ul  
 #C41243-AF680 100ul #C41243-AF750 100ul #C41243-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
 Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |  |
|-----------------------|--|
| Product Name          | NK-1R Polyclonal Conjugated Antibody   |
| Host Species          | Rabbit   |
| Clonality             | Polyclonal   |
| Specificity           | NK-1R Polyclonal Antibody detects endogenous levels of NK-1R protein.  |
| Immunogen Description | Synthesized peptide derived from the Internal region of human NK-1R.   |
| Conjugates            | Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750   |
| Other Names           | TACR1; NK1R; TAC1R; Substance-P receptor; SPR; NK-1 receptor; NK-1R; Tachykinin receptor 1   |
| Accession No.         | Swiss-Prot#:P25103NCBI Gene ID:6869  |
| Uniprot               | P25103   |
| GeneID                | 6869;  |
| Excitation Emission   | AF350: 346nm/442nm<br>AF405: 401nm/421nm<br>AF488: 493nm/519nm<br>AF555: 555nm/565nm<br>AF594: 591nm/614nm<br>AF647: 651nm/667nm<br>AF680: 679nm/702nm<br>AF750: 749nm/775nm |
| Calculated MW         | 50   |
| Concentration         | 0.5 mg/ml  |
| Formulation           | 0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide  |
| Storage               | Store at 4°C in dark for 6 months  |

## Application Details

## Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

---

Note: This product is for in vitro research use only