

Olfactory receptor 4D1 Polyclonal Conjugated Antibody

Catalog No: #C41272

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #C41272-AF350 100ul #C41272-AF405 100ul #C41272-AF488 100ul

#C41272-AF555 100ul #C41272-AF594 100ul #C41272-AF647 100ul

#C41272-AF680 100ul #C41272-AF750 100ul #C41272-Biotin 100ul

Description

Product Name	Olfactory receptor 4D1 Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Specificity	Olfactory receptor 4D1 Polyclonal Antibody detects endogenous levels of Olfactory receptor 4D1 protein.
Immunogen Description	Synthesized peptide derived from the C-terminal region of human Olfactory receptor 4D1.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	OR4D1; OR4D3; Olfactory receptor 4D1; Olfactory receptor 4D3; Olfactory receptor TPCR16
Accession No.	Swiss-Prot#:Q15615NCBI Gene ID:26689
Uniprot	Q15615
GeneID	26689;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	34
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only