

# Olfactory receptor 4K14 Polyclonal Conjugated Antibody

Catalog No: #C41273

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #C41273-AF350 100ul #C41273-AF405 100ul #C41273-AF488 100ul

#C41273-AF555 100ul #C41273-AF594 100ul #C41273-AF647 100ul

#C41273-AF680 100ul #C41273-AF750 100ul #C41273-Biotin 100ul

## Description

Product Name	Olfactory receptor 4K14 Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Specificity	Olfactory receptor 4K14 Polyclonal Antibody detects endogenous levels of Olfactory receptor 4K14 protein.
Immunogen Description	Synthesized peptide derived from the C-terminal region of human Olfactory receptor 4K14.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	OR4K14; Olfactory receptor 4K14; Olfactory receptor OR14-22
Accession No.	Swiss-Prot#:Q8NGD5NCBI Gene ID:122740
Uniprot	Q8NGD5
GeneID	122740;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	35
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

## Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

---

Note: This product is for in vitro research use only