

PKC Polyclonal Conjugated Antibody

Catalog No: #C41347



Package Size: #C41347-AF350 100ul #C41347-AF405 100ul #C41347-AF488 100ul
 #C41347-AF555 100ul #C41347-AF594 100ul #C41347-AF647 100ul
 #C41347-AF680 100ul #C41347-AF750 100ul #C41347-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	PKC Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Specificity	PKC Polyclonal Antibody detects endogenous levels of PKC protein.
Immunogen Description	Synthesized peptide derived from human PKC around the non-phosphorylation site of T497.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	PRKCA; PKCA; PRKACA; Protein kinase C alpha type; PKC-A; PKC-alpha; PRKCB; PKCB; PRKCB1; Protein kinase C beta type; PKC-B; PKC-beta; PRKCD; Protein kinase C delta type; Tyrosine-protein kinase PRKCD; nPKC-delta; PRKCE; PKCE; Protein kinase
Accession No.	Swiss-Prot#:P17252/P05771/Q05655/Q02156/P05129/P24723/Q04759/Q05513NCBI Gene ID:5578/5579/5580/5581/5582/5583/5588/5590
Uniprot	P17252
GeneID	5578;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	83
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only