

POTE-14/22 Polyclonal Conjugated Antibody

Catalog No: #C41356



Package Size: #C41356-AF350 100ul #C41356-AF405 100ul #C41356-AF488 100ul
 #C41356-AF555 100ul #C41356-AF594 100ul #C41356-AF647 100ul
 #C41356-AF680 100ul #C41356-AF750 100ul #C41356-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	POTE-14/22 Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Specificity	POTE-14/22 Polyclonal Antibody detects endogenous levels of POTE-14/22 protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human POTE-14/22.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	POTEH; A26C3; ACTBL1; POTE22; POTE ankyrin domain family member H; ANKRD26-like family C member 3; Prostate; ovary, testis-expressed protein on chromosome 22; POTE-22; POTE G; A26C2; POTE14; POTE ankyrin domain family member G; ANKRD26-like
Accession No.	Swiss-Prot#:Q6S545/Q6S5H5NCBI Gene ID:23784/404785
Uniprot	Q6S545
GeneID	23784;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	61
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250
 AF405 conjugated: most applications: 1: 50 - 1: 250
 AF488 conjugated: most applications: 1: 50 - 1: 250
 AF555 conjugated: most applications: 1: 50 - 1: 250
 AF594 conjugated: most applications: 1: 50 - 1: 250
 AF647 conjugated: most applications: 1: 50 - 1: 250
 AF680 conjugated: most applications: 1: 50 - 1: 250
 AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only