

Ribosomal Protein L10L Polyclonal Conjugated Antibody

Catalog No: #C41399

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #C41399-AF350 100ul #C41399-AF405 100ul #C41399-AF488 100ul

#C41399-AF555 100ul #C41399-AF594 100ul #C41399-AF647 100ul

#C41399-AF680 100ul #C41399-AF750 100ul #C41399-Biotin 100ul

Description

Product Name	Ribosomal Protein L10L Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Specificity	Ribosomal Protein L10L Polyclonal Antibody detects endogenous levels of Ribosomal Protein L10L protein.
Immunogen Description	Synthesized peptide derived from the C-terminal region of human Ribosomal Protein L10L.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	RPL10L; 60S ribosomal protein L10-like
Accession No.	Swiss-Prot#:Q96L21NCBI Gene ID:140801
Uniprot	Q96L21
GeneID	140801;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	20
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only