

Sg III Polyclonal Conjugated Antibody

Catalog No: #C41557



Package Size: #C41557-AF350 100ul #C41557-AF405 100ul #C41557-AF488 100ul
 #C41557-AF555 100ul #C41557-AF594 100ul #C41557-AF647 100ul
 #C41557-AF680 100ul #C41557-AF750 100ul #C41557-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	Sg III Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Specificity	Sg III Polyclonal Antibody detects endogenous levels of Sg III protein.
Immunogen Description	Synthesized peptide derived from the N-terminal region of human Sg III.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	SCG3; Secretogranin-3; Secretogranin III; SgIII
Accession No.	Swiss-Prot#:Q8WXD2NCBI Gene ID:29106
Uniprot	Q8WXD2
GeneID	29106;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	55
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only