

IL-5R α Polyclonal Conjugated Antibody

Catalog No: #C41622



Package Size: #C41622-AF350 100ul #C41622-AF405 100ul #C41622-AF488 100ul
 #C41622-AF555 100ul #C41622-AF594 100ul #C41622-AF647 100ul
 #C41622-AF680 100ul #C41622-AF750 100ul #C41622-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	IL-5R α Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Specificity	IL-5R α Polyclonal Antibody detects endogenous levels of IL-5R α protein.
Immunogen Description	Synthesized peptide derived from the N-terminal region of human IL-5R α .
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	IL5RA; IL5R; Interleukin-5 receptor subunit alpha; IL-5 receptor subunit alpha; IL-5R subunit alpha; IL-5R-alpha; IL-5RA; CDw125; CD125
Accession No.	Swiss-Prot#:Q01344NCBI Gene ID:3568
Uniprot	Q01344
GeneID	3568;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	47
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only