

IL-1RI Polyclonal Conjugated Antibody

Catalog No: #C41670



Package Size: #C41670-AF350 100ul #C41670-AF405 100ul #C41670-AF488 100ul
 #C41670-AF555 100ul #C41670-AF594 100ul #C41670-AF647 100ul
 #C41670-AF680 100ul #C41670-AF750 100ul #C41670-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	IL-1RI Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Specificity	IL-1RI Polyclonal Antibody detects endogenous levels of IL-1RI protein.
Immunogen Description	Synthesized peptide derived from the N-terminal region of human IL-1RI.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	IL1R1; IL1R; IL1RA; IL1RT1; Interleukin-1 receptor type 1; IL-1R-1; IL-1RT-1; IL-1RT1; CD121 antigen-like family member A; Interleukin-1 receptor alpha; IL-1R-alpha; Interleukin-1 receptor type I; p80; CD121a
Accession No.	Swiss-Prot#:P14778NCBI Gene ID:3554
Uniprot	P14778
GeneID	3554;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	68
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only