

HNF-4 α/γ Polyclonal Conjugated Antibody

Catalog No: #C41720



Package Size: #C41720-AF350 100ul #C41720-AF405 100ul #C41720-AF488 100ul
 #C41720-AF555 100ul #C41720-AF594 100ul #C41720-AF647 100ul
 #C41720-AF680 100ul #C41720-AF750 100ul #C41720-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	HNF-4 α/γ Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Specificity	HNF-4 α/γ Polyclonal Antibody detects endogenous levels of HNF-4 α/γ protein.
Immunogen Description	Synthesized peptide derived from human HNF-4 α/γ around the non-acetylation site of K127.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	HNF4A; HNF4; NR2A1; TCF14; Hepatocyte nuclear factor 4-alpha; HNF-4-alpha; Nuclear receptor subfamily 2 group A member 1; Transcription factor 14; TCF-14;
Accession No.	Swiss-Prot#:P41235NCBI Gene ID:3172
Uniprot	P41235
GeneID	3172;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	52
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only