

## CacyBP Polyclonal Conjugated Antibody

Catalog No: #C41775



Package Size: #C41775-AF350 100ul #C41775-AF405 100ul #C41775-AF488 100ul  
 #C41775-AF555 100ul #C41775-AF594 100ul #C41775-AF647 100ul  
 #C41775-AF680 100ul #C41775-AF750 100ul #C41775-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
 Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	CacyBP Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Specificity	CacyBP Polyclonal Antibody detects endogenous levels of CacyBP protein.
Immunogen Description	Synthesized peptide derived from the C-terminal region of human CacyBP.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	CACYBP; S100A6BP; SIP; PNAS-107; Calcyclin-binding protein; CacyBP; hCacyBP; S100A6-binding protein; Siah-interacting protein
Accession No.	Swiss-Prot#:Q9HB71NCBI Gene ID:27101
Uniprot	Q9HB71
GeneID	27101;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	25
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

## Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

---

Note: This product is for in vitro research use only