

BMPR1A Polyclonal Conjugated Antibody

Catalog No: #C41897



Package Size: #C41897-AF350 100ul #C41897-AF405 100ul #C41897-AF488 100ul
 #C41897-AF555 100ul #C41897-AF594 100ul #C41897-AF647 100ul
 #C41897-AF680 100ul #C41897-AF750 100ul #C41897-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	BMPR1A Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Specificity	BMPR1A Polyclonal Antibody detects endogenous levels of BMPR1A protein.
Immunogen Description	Synthesized peptide derived from the N-terminal region of human BMPR1A.
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	BMPR1A; ACVRLK3; ALK3; Bone morphogenetic protein receptor type-1A; BMP type-1A receptor; BMPR-1A; Activin receptor-like kinase 3; ALK-3; Serine/threonine-protein kinase receptor R5; SKR5; CD292
Accession No.	Swiss-Prot#:P36894NCBI Gene ID:657
Uniprot	P36894
GeneID	657;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	60
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Note: This product is for in vitro research use only