

## Interferon gamma Conjugated Monoclonal Antibody

Catalog No: #C42034



Package Size: #C42034-AF350 100ul #C42034-AF405 100ul #C42034-AF488 100ul  
 #C42034-AF555 100ul #C42034-AF594 100ul #C42034-AF647 100ul  
 #C42034-AF680 100ul #C42034-AF750 100ul #C42034-Biotin 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
 Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |  |
|-----------------------|--|
| Product Name          | Interferon gamma Conjugated Monoclonal Antibody  |
| Host Species          | Mouse  |
| Clonality             | Monoclonal   |
| Species Reactivity    | Hu   |
| Specificity           | specific for Human Interferon gamma denatured forms  |
| Immunogen Description | Recombinant Human Interferon gamma protein   |
| Conjugates            | Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750   |
| Other Names           | IFNG,Immune interferon,IFN-gamma   |
| Accession No.         | Swiss-Prot#:P01579   |
| Uniprot               | P01579   |
| GeneID                | 3458;  |
| Excitation Emission   | AF350: 346nm/442nm<br>AF405: 401nm/421nm<br>AF488: 493nm/519nm<br>AF555: 555nm/565nm<br>AF594: 591nm/614nm<br>AF647: 651nm/667nm<br>AF680: 679nm/702nm<br>AF750: 749nm/775nm |
| Formulation           | 0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide  |
| Storage               | Store at 4°C in dark for 6 months  |

## Application Details

## Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

## Background

---

Produced by lymphocytes activated by specific antigens or mitogens. IFN-gamma, in addition to having antiviral activity, has important immunoregulatory functions. It is a potent activator of macrophages, it has antiproliferative effects on transformed cells and it can potentiate the antiviral and antitumor effects of the type I interferons.

---

Note: This product is for in vitro research use only