

Kai1 antigen Polyclonal Conjugated Antibody

Catalog No: #C42112



Package Size: #C42112-AF350 100ul #C42112-AF405 100ul #C42112-AF488 100ul
 #C42112-AF555 100ul #C42112-AF594 100ul #C42112-AF647 100ul
 #C42112-AF680 100ul #C42112-AF750 100ul #C42112-Biotin 100ul

Orders: order@signalwayantibody.com
 Support: tech@signalwayantibody.com

Description

Product Name	Kai1 antigen Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Ms
Specificity	The antibody detects endogenous level of total Kai1 antigen polyclonal antibody.
Immunogen Description	Recombinant mouse Kai1 antigen protein
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	C33 antigen, IA4, Inducible membrane protein R2, Metastasis suppressor Kangai-1 homolog, CD82, Cd82, Kai1
Accession No.	Swiss-Prot#: P40237
Uniprot	P40237
GeneID	12521;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	30
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Background

Associates with CD4 or CD8 and delivers costimulatory signals for the TCR/CD3 pathway.

Note: This product is for in vitro research use only