

Cytochrome c oxidase subunit 4 isoform 1, mitochondrial Polyclonal Conjugated Antibody

Catalog No: #C42116

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #C42116-AF350 100ul #C42116-AF405 100ul #C42116-AF488 100ul

#C42116-AF555 100ul #C42116-AF594 100ul #C42116-AF647 100ul

#C42116-AF680 100ul #C42116-AF750 100ul #C42116-Biotin 100ul

Description

Product Name	Cytochrome c oxidase subunit 4 isoform 1, mitochondrial Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Hu
Specificity	The antibody detects endogenous level of total Cytochrome c oxidase subunit 4 isoform 1, mitochondrial polyclonal antibody.
Immunogen Description	Recombinant human Cytochrome c oxidase subunit 4 isoform 1, mitochondrial protein
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Cytochrome c oxidase polypeptide IV, Cytochrome c oxidase subunit IV isoform 1, COX IV-1, COX4I1, COX4
Accession No.	Swiss-Prot#:P13073
Uniprot	P13073
GeneID	1327;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	20
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Background

This protein is one of the nuclear-coded polypeptide chains of cytochrome c oxidase, the terminal oxidase in mitochondrial electron transport.

Note: This product is for in vitro research use only