

Protein-lysine 6-oxidas Polyclonal Conjugated Antibody

Catalog No: #C42240

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #C42240-AF350 100ul #C42240-AF405 100ul #C42240-AF488 100ul

#C42240-AF555 100ul #C42240-AF594 100ul #C42240-AF647 100ul

#C42240-AF680 100ul #C42240-AF750 100ul #C42240-Biotin 100ul

Description

Product Name	Protein-lysine 6-oxidas Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Rt
Specificity	The antibody detects endogenous level of total Protein-lysine 6-oxidas polyclonal antibody.
Immunogen Description	Recombinant Rat Protein-lysine 6-oxidas protein
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Lysyl oxidase Lox
Accession No.	Swiss-Prot#:P28301
Uniprot	P28301
GeneID	16948;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	47
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Background

Responsible for the post-translational oxidative deamination of peptidyl lysine residues in precursors to fibrous collagen and elastin.

Note: This product is for in vitro research use only