

Non-secretory ribonuclease Polyclonal Conjugated Antibody

Catalog No: #C42307

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #C42307-AF350 100ul #C42307-AF405 100ul #C42307-AF488 100ul

#C42307-AF555 100ul #C42307-AF594 100ul #C42307-AF647 100ul

#C42307-AF680 100ul #C42307-AF750 100ul #C42307-Biotin 100ul

Description

Product Name	Non-secretory ribonuclease Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Hu
Specificity	The antibody detects endogenous level of total Non-secretory ribonuclease polyclonal antibody.
Immunogen Description	Recombinant human Non-secretory ribonuclease protein
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Eosinophil-derived neurotoxin,RNase Upl-2,Ribonuclease 2,Ribonuclease US,RNASE2,EDN, RNS2
Accession No.	Swiss-Prot#:P10153
Uniprot	P10153
GeneID	6036;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Background

This is a non-secretory ribonuclease. It is a pyrimidine specific nuclease with a slight preference for U. Cytotoxic and helminthotoxin. Selectively chemotactic for dendritic cells. Possesses a wide variety of biological activities.

Note: This product is for in vitro research use only