

KITH_HHV11 Antibody

Catalog No: #45481



Package Size: #45481-1 50ul #45481-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	KITH_HHV11 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	KITH_HHV11 Antibody detects endogenous levels of total KITH_HHV11
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human KITH_HHV11
Target Name	KITH_HHV11
Other Names	Thymidine kinase
Accession No.	Swiss-Prot#: P0DTH5/Q9QNF7NCBI Gene ID:24271467
Uniprot	P0DTH5/Q9QNF7
GeneID	24271467;
Calculated MW	71 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Note: This product is for in vitro research use only