

UV excision repair protein RAD23 homolog A Polyclonal Conjugated Antibody

Catalog No: #C42402

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #C42402-AF350 100ul #C42402-AF405 100ul #C42402-AF488 100ul

#C42402-AF555 100ul #C42402-AF594 100ul #C42402-AF647 100ul

#C42402-AF680 100ul #C42402-AF750 100ul #C42402-Biotin 100ul

Description

Product Name	UV excision repair protein RAD23 homolog A Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Hu
Specificity	The antibody detects endogenous level of total UV excision repair protein RAD23 homolog A polyclonal antibody.
Immunogen Description	Recombinant human UV excision repair protein RAD23 homolog A protein
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	HR23A hHR23A
Accession No.	Swiss-Prot#:P54725
Uniprot	P54725
GeneID	5886;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	40
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Background

Multiubiquitin chain receptor involved in modulation of proteasomal degradation. Binds to 'Lys-48'-linked polyubiquitin chains in a length-dependent manner and with a lower affinity to 'Lys-63'-linked polyubiquitin chains. Proposed to be capable to bind simultaneously to the 26S proteasome and to polyubiquitinated substrates and to deliver ubiqui

Note: This product is for in vitro research use only