

C-C motif chemokine 24 Polyclonal Conjugated Antibody

Catalog No: #C42444

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #C42444-AF350 100ul #C42444-AF405 100ul #C42444-AF488 100ul

#C42444-AF555 100ul #C42444-AF594 100ul #C42444-AF647 100ul

#C42444-AF680 100ul #C42444-AF750 100ul #C42444-Biotin 100ul

Description

Product Name	C-C motif chemokine 24 Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Hu
Specificity	The antibody detects endogenous level of total C-C motif chemokine 24 polyclonal antibody.
Immunogen Description	Recombinant human C-C motif chemokine 24 protein
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	MPIF2, SCYA24,CCL24,CK-beta-6,Eosinophil chemotactic protein 2,Eotaxin-2,Myeloid progenitor inhibitory factor 2,MPIF-2,Small-inducible cytokine A24
Accession No.	Swiss-Prot#:O00175
Uniprot	O00175
GeneID	6369;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Background

Chemotactic for resting T-lymphocytes, and eosinophils. Has lower chemotactic activity for neutrophils but none for monocytes and activated lymphocytes. Is a strong suppressor of colony formation by a multipotential hematopoietic progenitor cell line. Binds to CCR3.

Note: This product is for in vitro research use only