

## CSRNP1 Polyclonal Conjugated Antibody

Catalog No: #C42663



Package Size: #C42663-AF350 100ul #C42663-AF405 100ul #C42663-AF488 100ul

#C42663-AF555 100ul #C42663-AF594 100ul #C42663-AF647 100ul

#C42663-AF680 100ul #C42663-AF750 100ul #C42663-Biotin 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

|                       |                                                                                                                                                                              |
|-----------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | CSRNP1 Polyclonal Conjugated Antibody                                                                                                                                        |
| Host Species          | Rabbit                                                                                                                                                                       |
| Clonality             | Polyclonal                                                                                                                                                                   |
| Species Reactivity    | Hu                                                                                                                                                                           |
| Specificity           | The antibody detects endogenous level of total CSRNP1 polyclonal antibody.                                                                                                   |
| Immunogen Description | Recombinant human Cysteine/serine-rich nuclear protein 1 protein(1-200aa)                                                                                                    |
| Conjugates            | Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750                                                                                                                       |
| Other Names           | CSRNP-1,Axin-1 up-regulated gene 1 protein,Protein URAX1,TGF-beta-induced apoptosis protein 3,TAIP-3,CSRNP1,AXUD1,TAIP3                                                      |
| Accession No.         | Swiss-Prot#:Q96S65                                                                                                                                                           |
| Uniprot               | Q96S65                                                                                                                                                                       |
| GeneID                | 64651;                                                                                                                                                                       |
| Excitation Emission   | AF350: 346nm/442nm<br>AF405: 401nm/421nm<br>AF488: 493nm/519nm<br>AF555: 555nm/565nm<br>AF594: 591nm/614nm<br>AF647: 651nm/667nm<br>AF680: 679nm/702nm<br>AF750: 749nm/775nm |
| Formulation           | 0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide                                                                                  |
| Storage               | Store at 4°C in dark for 6 months                                                                                                                                            |

## Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

## Background

---

Binds to the consensus sequence 5'-AGAGTG-3' and has transcriptional activator activity (By similarity). May have a tumor-suppressor function. May play a role in apoptosis.

---

Note: This product is for in vitro research use only