

SPT4H Antibody

Catalog No: #45510

Package Size: #45510-1 50ul #45510-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	SPT4H Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	SPT4H Antibody detects endogenous levels of total SPT4H
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human SPT4H
Target Name	SPT4H
Other Names	SUPT4H1, DSIF p14, SPT4, SPT4H, SUPT4H, DSIF small subunit, HSPT4
Accession No.	Swiss-Prot#: P63272NCBI Gene ID:6827
Uniprot	P63272
GeneID	6827;
Calculated MW	13 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Note: This product is for in vitro research use only