

Ubiquitin-conjugating enzyme E2 Q2 Polyclonal Conjugated Antibody

Catalog No: #C42667

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #C42667-AF350 100ul #C42667-AF405 100ul #C42667-AF488 100ul

#C42667-AF555 100ul #C42667-AF594 100ul #C42667-AF647 100ul

#C42667-AF680 100ul #C42667-AF750 100ul #C42667-Biotin 100ul

Description

Product Name	Ubiquitin-conjugating enzyme E2 Q2 Polyclonal Conjugated Antibody
Host Species	Rabbit
Clonality	Polyclonal
Species Reactivity	Hu
Specificity	The antibody detects endogenous level of total Ubiquitin-conjugating enzyme E2 Q2 polyclonal antibody.
Immunogen Description	Recombinant human Ubiquitin-conjugating enzyme E2 Q2 protein
Conjugates	Biotin AF350 AF405 AF488 AF555 AF594 AF647 AF680 AF750
Other Names	Ubiquitin carrier protein Q2 Ubiquitin-protein ligase Q2 Gene names UBE2Q2
Accession No.	Swiss-Prot#:Q8WVN8
Uniprot	Q8WVN8
GeneID	92912;
Excitation Emission	AF350: 346nm/442nm AF405: 401nm/421nm AF488: 493nm/519nm AF555: 555nm/565nm AF594: 591nm/614nm AF647: 651nm/667nm AF680: 679nm/702nm AF750: 749nm/775nm
Calculated MW	43
Formulation	0.01M Sodium Phosphate, 0.25M NaCl, pH 7.6, 5mg/ml Bovine Serum Albumin, 0.02% Sodium Azide
Storage	Store at 4°C in dark for 6 months

Application Details

Suggested Dilution:

AF350 conjugated: most applications: 1: 50 - 1: 250

AF405 conjugated: most applications: 1: 50 - 1: 250

AF488 conjugated: most applications: 1: 50 - 1: 250

AF555 conjugated: most applications: 1: 50 - 1: 250

AF594 conjugated: most applications: 1: 50 - 1: 250

AF647 conjugated: most applications: 1: 50 - 1: 250

AF680 conjugated: most applications: 1: 50 - 1: 250

AF750 conjugated: most applications: 1: 50 - 1: 250

Biotin conjugated: working with enzyme-conjugated streptavidin, most applications: 1: 50 - 1: 1,000

Background

Accepts ubiquitin from the E1 complex and catalyzes its covalent attachment to other proteins. In vitro catalyzes 'Lys-48'-linked polyubiquitination.

Note: This product is for in vitro research use only