

KKLC1 Antibody

Catalog No: #45529

Package Size: #45529-1 50ul #45529-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	KKLC1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	KKLC1 Antibody detects endogenous levels of total KKLC1
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human KKLC1
Target Name	KKLC1
Other Names	CXorf61, KKLC1
Accession No.	Swiss-Prot#: Q5H943NCBI Gene ID:203413
Uniprot	Q5H943
GeneID	203413;
Calculated MW	13 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Note: This product is for in vitro research use only