

TNFAIP8L1 Antibody

Catalog No: #45539

Package Size: #45539-1 50ul #45539-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	TNFAIP8L1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	TNFAIP8L1 Antibody detects endogenous levels of total TNFAIP8L1
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human TNFAIP8L1
Target Name	TNFAIP8L1
Other Names	Oxidative stress-regulated gene-beta
Accession No.	Swiss-Prot#: Q8WVP5NCBI Gene ID:126282
Uniprot	Q8WVP5
GeneID	126282;
Calculated MW	21 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Note: This product is for in vitro research use only