

SRRM2 Antibody

Catalog No: #45566

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SRRM2 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	SRRM2 Antibody detects endogenous levels of total SRRM2
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human SRRM2
Target Name	SRRM2
Other Names	SRRM2, 300 kDa nuclear matrix antigen, 300-KD, Cwc21, KIAA0324, RNA binding protein, SRm300, TaxREB803, SRL300, CWF21
Accession No.	Swiss-Prot#: Q9UQ35NCBI Gene ID:23524
Uniprot	Q9UQ35
GeneID	23524;
Calculated MW	300kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Note: This product is for in vitro research use only