

THRAP3 Antibody

Catalog No: #45570

Package Size: #45570-1 50ul #45570-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	THRAP3 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	THRAP3 Antibody detects endogenous levels of total THRAP3
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human THRAP3
Target Name	THRAP3
Other Names	THRAP3, TRAP150
Accession No.	Swiss-Prot#: Q9Y2W1NCBI Gene ID:9967
Uniprot	Q9Y2W1
GeneID	9967;
Calculated MW	109kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western Blot analysis of HeLa lysates, using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only