

CD320 Antibody

Catalog No: #45581

Package Size: #45581-1 50ul #45581-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	CD320 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	CD320 Antibody detects endogenous levels of total CD320
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human CD320
Target Name	CD320
Other Names	CD320, 8D6 antigen, 8D6, 8D6A, CD320 antigen, CD320 molecule, FDC-SM-8D6, FDC-signaling molecule 8D6, Transcobalamin receptor, TCBLR
Accession No.	Swiss-Prot#: Q9NPF0NCBI Gene ID:51293
Uniprot	Q9NPF0
GeneID	51293;
Calculated MW	29kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis of lysates from SW480 cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only