

## Cytochrome c Oxidase 1 Antibody

Catalog No: #45586



Package Size: #45586-1 50ul #45586-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Cytochrome c Oxidase 1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu Ms
Specificity	Cytochrome c Oxidase 1 Antibody detects endogenous levels of total Cytochrome c Oxidase 1
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human Cytochrome c Oxidase 1
Target Name	Cytochrome c Oxidase 1
Other Names	MT-COI, COX1, COI, MT-C01, MTCO1
Accession No.	Swiss-Prot#: P00395NCBI Gene ID:4512
Uniprot	P00395
GeneID	4512;
Calculated MW	57 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000-3000

Note: This product is for in vitro research use only