

Cytochrome P450 27C1 Antibody

Catalog No: #45588

Package Size: #45588-1 50ul #45588-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Cytochrome P450 27C1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	Cytochrome P450 27C1 Antibody detects endogenous levels of total Cytochrome P450 27C1
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human Cytochrome P450 27C1
Target Name	Cytochrome P450 27C1
Other Names	CYP27C1, Cytochrome P450 27C1
Accession No.	Swiss-Prot#: Q4G0S4NCBI Gene ID:339761
Uniprot	Q4G0S4
GeneID	339761;
Calculated MW	43 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only