

Cytochrome P450 4F22 Antibody

Catalog No: #45595

Package Size: #45595-1 50ul #45595-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	Cytochrome P450 4F22 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Immunogen affinity purified
Applications	WB
Species Reactivity	Hu
Specificity	Cytochrome P450 4F22 Antibody detects endogenous levels of total Cytochrome P450 4F22
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human Cytochrome P450 4F22
Target Name	Cytochrome P450 4F22
Other Names	CYP4F22, Cytochrome P450 4F22, LI3
Accession No.	Swiss-Prot#: Q6NT55NCBI Gene ID:126410
Uniprot	Q6NT55
GeneID	126410;
Calculated MW	62 kDa
Concentration	1.0mg/mL
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western Blot analysis of HEK293 lysates, using primary antibody at 1:1000 dilution.

Note: This product is for in vitro research use only